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
FOR LÆREREN:

EQUIPMENT TIPS

1. Setting an appropriate release pressure for the Projectile Launcher is critical. You will need to give students additional details to help them set it. Turn the Range knob clockwise for higher pressure and higher muzzle speed and counter-clockwise for lower pressure and lower muzzle speed.

a. To set the pressure, we recommend turning the Range knob to the highest pressure (all the way to the right) and then turn it counter-clockwise a couple of times.

b. Build pressure using the pump until the pressure stabilizes. To do this, when you pump the hand pump, you will hear a small release sound when the pressure is reached. Keep Physics with Vernier © Vernier Software & Technology 8B - 1 I INSTRUCTOR INFORMATION pumping until you hear at least three small release sounds and then wait for five seconds so that the pressure is fully stabilized.

c. Look at the valve on the pump. Pressures between 40 and 50 psi work well for this experiment. If the pressure readings are above or below this range, adjust the Range knob and build pressure again.

d. You may want to set the pressure for an appropriate firing distance before your students see the apparatus and instruct students not to adjust the Range knob.

e. Impact point predictions will be changed if the Projectile Launcher moves during trials. Consider clamping the Projectile Launcher to the table.

f. When using the Time of Flight Pad (order code: TOF-VPL) with the Projectile Launcher, be sure to select Time of Flight Pad during data collection setup in LabQuest and Graphical Analysis.

DATA COLLECTION AND ANALYSIS TIPS

1. Encourage students to algebraically solve the projectile motion problem in the Analysis section (that is, without numbers), since they will evaluate the function three times.

2. There are a number of ways to account for variation in the speed measurement in the range prediction, Δx . Depending on the mathematical sophistication of your students, use either the average \pm standard deviation or the minimum and maximum speeds as extreme values for the speed. Calculate the range for each extreme value, creating a range of predictions.

3. If you hold a competition among your students in predicting the impact point, note that fine distinctions are often statistically impossible to defend. For example, suppose one group predicts a range of 0.78 ± 0.03 m, and another 0.72 ± 0.04 m. The actual impact points are 0.76 and 0.71 m, respectively. It would be wrong to say that the second group was “closer” in their experiment, since the impact was within the predicted range for both groups. It was simple random variation that caused the second group to have a strike that is closer to the center of the predicted range than

the first. If the trials were repeated, it is possible that the first group's deviation from the center of the impact range would be closer. Both groups are equally successful, and to rank the groups as second and first cannot be defended statistically. Give the groups equal prizes or bonus points.

ANSWERS TO PRELIMINARY QUESTIONS

1. Some students will incorrectly predict that the penny that is moving faster (horizontally) will land after (or before) the falling penny. Since the horizontal and vertical components of the motion are independent, they will land at the same time.
2. Students should note that the sounds of each penny striking the floor are very close together. Some students will hold on to their preconception, insisting that one landed before the other, even by a tiny bit. Have students discuss what evidence will sufficiently convince them otherwise. 8B - 2 I Physics with Vernier Experiment 8B Projectile Motion (Projectile Launcher) Physics with Vernier 8B - 3 I
3. There will be some variability in the landings. One should not consistently land before the other.
4. If students think the horizontal speed makes a difference, then the increased horizontal speed should make more of a difference.
5. If there is a difference, the height change should exaggerate it. One should land twice as much before as the other.
6. No, the horizontal speed does not matter because the pennies seem to hit the floor simultaneously.
7. Since the pennies hit the floor together, they take the same time to fall

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ANSWERS TO ANALYSIS QUESTIONS

1. Any calculation using uncertain measurements is itself uncertain, so the prediction cannot be exact. A range of predicted values is more appropriate than a single value.
2. Explain to your students that they should predict a target range rather than a target spot. The range depends on the uncertainty in measurements used in the calculation. As long as the actual impact point is somewhere between the minimum and maximum distances, where the ball strikes within the range depends on random variations.
3. In addition to variations in the speed measurement (which is itself affected by the Photogate spacing), the table height, h , may be different from the measured value.
4. No. No account of air resistance is made in this analysis. If air resistance is important then the ball falls short of the predicted impact point.